

NALGENE® 2250 Series

Autoclavable Polypropylene Carboys with Handles



Popular round shape with handles and 83B closure, autoclavable Polypropylene.

Excellent for large batches of culture media, distilled water and other solutions. Graduated in 1-gallon and 5-liter increments. Accepts No. 13-1/2 rubber stopper. NOTE: For best results, autoclave using a properly vented closure.

Not recommended for use with hazardous materials.

BIOPHARMACEUTICAL PRODUCTS

Features/Benefits:

- **Autoclavable Polypropylene**
- **Unbreakable, chemically inert and ultra-pure**
- **Non-stick; easy to clean**
- **Wide shoulder handles for easy carrying and pouring**
- **One-Piece construction minimizes risk of leakage**
- **Graduated in gallons and liters**
- **83B closure and neck insures a tight closure**

Cat. No.	2250-0020	2250-0050	2250-0130
Nom. Should. Cap., L	10	20	50
Nom. Should. Cap., gal.	2-1/2	5-1/2	13
Brim Cap., L (approx.)	12.5	24	55
No. per Case	6	4	1

See next page for properties & characteristics Polypropylene

sani-tech
west

321 Irving Drive
Oxnard, CA 93030
(800) 726-4835
www.sani-techwest.com

Page 1

Physical Properties of Polypropylene

Material	PP
Max Use Temperature	135°C
Heat Deflection Temperature ¹	107°C
Brittleness Temperature ¹¹	0°C
Transparency	Translucent
Microwavability ²	Yes
Autoclaving ⁴	Yes
Gas ⁴	Yes
Dry Heat ⁴	No
Radiation ⁴	No
Disinfectants	Yes
Specific Gravity	0.9
Flexibility	rigid
Permeability N2	48 cc•mil/100in ² •24hr•atm
Permeability O2	240 cc•mil/100in ² •24hr•atm
Permeability CO2	800 cc•mil/100in ² •24hr•atm
Permeability N2 (metric)	744 cc•mm/m ² •24hr•Bar
Permeability O2 (metric)	3,720 cc•mm/m ² •24hr•Bar
Permeability CO2 (metric)	12,400 cc•mm/m ² •24hr•Bar
WVTR	3.9 g•mm/m ² •24hr•Bar @ 37°C and 90% humidity
Water Adsorption	<0.02
Non-cytotoxicity ⁶	Yes
Food and Bev use Rating ⁷	Yes
CFR21	177.152

- Heat Deflection Temperature is the temperature at which a bar deflects 0.01" at 66 psig (ASTM D648). Materials may be used above Heat Deflection temperatures in non-stress applications; see Max. Use Temp.
- Ratings based on 5-minute tests using 600 watts of power on exposed, empty labware. CAUTION: Do not exceed Max. Use Temp., or expose labware to chemicals which heating cause to attack the plastic or be rapidly absorbed.
- Plastic will absorb heat.
- STERILIZATION:
 - Autoclaving (121°C, 15 psig for 20 minutes) -- Clean and rinse items with distilled water before autoclaving. (Always completely disengage thread before autoclaving.) Certain chemicals which have no appreciable effect on resins at room temperature may cause deterioration at autoclaving temperatures unless removed with distilled water beforehand.
 - Gas -- Ethylene Oxide, formaldehyde, hydrogen peroxide
 - Dry Heat (160°C, 120 minutes)
 - Disinfectants -- Benzalkonium chloride, formalin/formaldehyde, ethanol, etc.
 - Radiation -- gamma irradiation at 25 kGy (2.5 MRad) with unstabilized plastic.
- Sterilizing reduces mechanical strength. Do not use PC vessels for vacuum applications if they have been autoclaved. Refer to Use and Care Guidelines for NALGENE Labware, for detailed information on sterilizing.
- "Yes" indicates the resin has been determined to be non-cytotoxic, based on USP and ASTM biocompatibility testing standards utilizing an MEM elution technique on a WI38 human diploid lung cell line.
- Resins meet requirements of CFR21 section of Food Additives Amendment of the Federal Food and Drug Act. End users are responsible for validation of compliance for specific containers used in conjunction with their particular packaging applications.
- Acceptable for aqueous foods only, at temperatures up to 121°C/250°F. Not sanctioned for use with alcoholic or fatty foods at any temperature.
- Acceptable for:
 - Nonacid, aqueous products; may contain salt, sugar or both (pH above 5.0)
 - Dairy products and modifications; oil-in-water emulsions, high or low fat
 - Moist bakery products with surface containing no free fat or oil
 - Dry solids with the surfaces containing no free fat or oil (no end-test required) and under all conditions as described in Table 2 of FDA Regulation 177.1520 except condition A - high temperature sterilization (e.g. over 100° C/212° F)
- Acceptable for:
 - Alcoholic foods containing not more than 15% (by volume) alcohol; fill and storage temperature not to exceed 49°C (120°F)
 - Non-alcoholic foods of hot fill to not exceed 82°C (180°F) and 49°C (120°F) in storage.
 - Not suitable for carbonated beverages or beer or packaging food requiring thermal processing.
 - Straight-sided jars, beakers and graduated cylinders only.
 - Acceptable for aqueous, oil, dairy, acidic, and alcoholic foods up to 71°C/160°F.
- The brittleness temperature is the temperature at which an item made from the resin may break or cracked if dropped. This is not the lowest use temperature if care is exercised in use and handling.

