

NEW

BIO-SIMPLEX™ POLYPROPYLENE CARBOYS



Saint-Gobain Performance Plastics introduces Bio-Simplex™ Polypropylene Carboys with EZ Top® Caps and single-use EZ Top® Assemblies for advanced liquid handling and solution storage for many pharma applications. Saint-Gobain offers four popular sizes of carboys: 5 L, 10 L, 20 L and 40 L manufactured from high quality polypropylene, which is resistant to chemicals and heating up to 275°F (135°C), and allows sterilization for use in pharmaceutical and clinical applications.

The carboy's body includes design features that allow the user to easily hold, lift, fill, and pour while the versatile EZ Top cap technology provides interchangeable and leak-proof top connections. The Bio-Simplex carboys are available with or without an autoclavable spigot. They also feature easy-to-read certified graduation marks with ± 5% accuracy, large molded marks for clearly identifiable volume and material. The Bio-Simplex system—carboy, cap, and spigot is autoclavable.

The Saint-Gobain Bio-Simplex Carboy Closure System is also offered with the EZ Top Closure Assembly for aseptic fluid transfer, sampling or sealing. The single-use EZ Top solutions are made of two different resins:

- ◆ C-Flex® thermoplastic elastomer
- ◆ Sani-Tech® high purity platinum-cured silicone

Features and Benefits:

- ◆ Slender ergonomic body design requiring less space in the laboratory
- ◆ Large top handles for full hand grip
- ◆ Bottom handles for secure grip when holding and transporting
- ◆ Minimum dead volume for maximum solution recovery
- ◆ Wide mouth for easy filling and cleaning
- ◆ Volume and material identification marks
- ◆ Large, easy to read graduation marks with ± 5% accuracy
- ◆ Open/closed version of EZ Top cap with replaceable gasket
- ◆ Perfect fit high-purity EZ Top Closure Assembly
- ◆ Autoclavable spigot with interchangeable option



Catalog No.	Description	Capacity L	Spigot	Cap Size mm	Qty/Pack
DCBo58o	Bio-Simplex Carboy with Closed EZ Top Cap	5	No	8o	1
DCBo58oS	Bio-Simplex Carboy with Closed EZ Top Cap	5	Yes	8o	1
DCB1o8o	Bio-Simplex Carboy with Closed EZ Top Cap	10	No	8o	1
DCB1o8oS	Bio-Simplex Carboy with Closed EZ Top Cap	10	Yes	8o	1
DCB2o12o	Bio-Simplex Carboy with Closed EZ Top Cap	20	No	12o	1
DCB2o12oS	Bio-Simplex Carboy with Closed EZ Top Cap	20	Yes	12o	1
DCB4o12o	Bio-Simplex Carboy with Closed EZ Top Cap	40	No	12o	1
DCB4o12oS	Bio-Simplex Carboy with Closed EZ Top Cap	40	Yes	12o	1

BIO-SIMPLEX™ CARBOY ACCESSORIES

Catalog No.	Description	Cap Size mm	Qty/Pack
DEZo8oCT	EZ Top Cap, 8o, Closed Top, White	8o	2
DEZo8oOT	EZ Top Cap, 8o, Open Top, White, with Closed Adapter	8o	2
DEZ12oCT	EZ Top Cap, 12o, Closed Top, White	12o	2
DEZ12oOT	EZ Top Cap, 12o, Open Top, White, with Closed Adapter	12o	2
DCBooSP	Spigot for Carboys		1
DEZ12oAC	Adapter Cone, EZ Top Cap, 12o, White		1

EZ TOP® CONTAINER CLOSURE ASSEMBLIES

The EZ Top Container Closure is used to create a complete container assembly when combined with cap and container. This economical assembly designed for single use can successfully replace stainless steel closures and eliminate costs such as cleaning, validation and assembly, and the contamination from repeated use.

C-FLEX® EZ TOP® CONTAINER CLOSURES

C-Flex EZ Top Container Closures are made from C-Flex pharmaceutical grade resin and tubing to ensure an unobstructed fluid path, extremely low level of extractables, and provide an elastomeric seal against glass, plastic, and metal surfaces. The one-piece construction eliminates product loss and contamination due to loose connection.

The C-Flex EZ Top Container Closure includes C-Flex tubing with standard plug and fitting, a Bio-Simplex™ Flex Joint, 0.2 µm vent filter, and cap.

Two or three port configurations are available to fit all sizes of Bio-Simplex carboys.

C-FLEX® EZ TOP® ASSEMBLIES FOR CARBOYS

Catalog No.	Container Size, L	Ports
CBO5L-2EZ	5	2
CB10L-2EZ	10	2
CB20L-2EZ	20	2
CBO5L-3EZ	5	3
CB10L-3EZ	10	3
CB20L-3EZ	20	3

Custom design of your C-Flex® EZ Top® container top assembly with factory installed 1-3 port options, C-Flex® tubing configurations, and sterilization are available upon request. Please contact customer service at 1-800-541-6880 or +1-727-531-4191 for assistance.



SANI-TECH® EZ TOP® CONTAINER CLOSURES

Sani-Tech EZ Top Container Closures are manufactured from fully characterized Class VI high-purity, platinum-cured silicone. Sani-Tech EZ Top Container Closure components provide an aseptic seal for handling various high purity fluids in biopharmaceutical applications. The smooth inner bore eliminates risk of product entrapment and buildup. Sani-Tech EZ Top Container Closures allow easy retrofit of your existing glass or plastic carboys for aseptic fluid transfer, sampling or sealing. Sani-Tech EZ Top Container Closures can be autoclaved and sterilized. It is non-cytotoxic and complies with USP Class VI standards.

- ◆ Silicone closure insert with or without silicone tubing or sanitary dip tubes
- ◆ Custom design capabilities for all components

Custom design of your Sani-Tech® EZ Top® container top assembly with factory installed 1-3 port options, Sani-Tech® Tubing configurations and sterilization are available upon request. Please contact customer service at 1-800-283-5963 or +1-518-283-5963 for assistance.

Saint-Gobain Performance Plastics-Clearwater

4451 110th Avenue North
Clearwater, FL 33762
Tel: (727) 531-4191
Toll-Free: (800) 541-6880
Fax: (727) 299-2569



C-Flex®, EZ Top® and Sani-Tech® are registered trademarks of Saint-Gobain Performance Plastics Corporation.



IMPORTANT: It is the user's responsibility to ensure the suitability and safety of Saint-Gobain Performance Plastics tubing for all intended uses. Laboratory and clinical tests must be conducted in accordance with applicable regulatory requirements in order to determine the safety and effectiveness for use of tubing in any particular application.

For a period of 6 months from the date of first sale, Saint-Gobain Performance Plastics Corporation warrants this product to be free from defects in materials and workmanship. Our only obligation will be to replace any portion proving defective, or at our option, to refund the purchase price thereof. User assumes all other risk, if any, including the risk of injury, loss or damage, direct or consequential, arising out of the use, misuse, or inability to use, this product. SAINT-GOBAIN PERFORMANCE PLASTICS CORPORATION SHALL NOT BE LIABLE FOR ANY SPECIAL, INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY KIND RESULTING FROM THE USE OF THIS PRODUCT. THIS WARRANTY IS IN LIEU OF THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR PARTICULAR PURPOSE, AND ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. No deviation is authorized.

Saint-Gobain Performance Plastics Corporation assumes no obligations or liability for any advice furnished by it, or for results obtained with respect to those products. All such advice is given and accepted at the buyer's risk.

Interphex-0311-CWTR

www.biopharm.saint-gobain.com

©2011 Saint-Gobain Performance Plastics Corporation



High Purity Liquid Handling Products
for Pharmaceutical & Biotech Industries
(800) 726-4835 | www.sani-techwest.com